



# A SAFE PAIR OF HANDS

University of Nottingham, April  
17th 2009

A SHIELD Scientific Education Programme



# Learning Objectives

- **Describe differences between gloves**
- **Distinguish three common reactions to wearing latex and synthetic gloves**
- **List the four components of gloves that can trigger reactions and sensitivities**



# Increased Reactions

- ① Dry chapped hands
- ① Cracking
- ① Eczema (reddening of the skin)
- ① Dermatitis (inflammation)



# New Symptoms

- ① Hives
- ① Rhinitis
- ① Asthma
- ① Conjunctivitis
- ① Anaphylaxis



# Glove Components

- ① Chemicals
- ① Proteins
- ① Powder
- ① Endotoxins



# Chemicals provide:

① Strength

② Stretchability

③ Storage preservation



# Latex Being Tapped from Rubber Tree – source of latex proteins







# Manufacturing Process - use of powder as a release agent







# Use of Powder as a donning agent - residues on hand





# Endotoxins (Pyrogens)

- Originate from raw latex materials or result from poor manufacturing practices
- Not destroyed by the sterilization process
- May be present on sterile surgical or cleanroom gloves



# Dermal Reactions

- Irritant Contact Dermatitis (also known as irritation, irritant dermatitis or dermatitis)
- Allergic Contact Dermatitis (also known as Type IV, Delayed Hypersensitivity or Chemical Allergy)
- Natural Rubber Latex Allergy (also known as Immediate Type Hypersensitivity, Protein Allergy or Type I)



# Irritant contact dermatitis

- Non-allergic condition
- Can effect everyone (i.e. latex and synthetic glove wearers)
- Both non-glove-associated and glove-associated irritant contact dermatitis



# Glove Related Irritant Contact Dermatitis

- ① Chemicals
- ① Powder
- ① Endotoxins
- ① Friction
- ① Air occlusion



# Non-Glove Related Irritant Contact Dermatitis

- ① Soaps
- ① Detergents
- ① Disinfectants
- ① Degreasing agents
- ① Ethylene Oxide
- ① Alcohol





# Acute Irritant Contact Dermatitis





# Irritant Contact Dermatitis – the Symptoms

- Itching
- Redness
- Inflammation
- Scaly appearance
- Burning sensation
- Thickened skin
- Blisters
- Small hard bumps
- Cracks
- Line of demarcation with glove induced irritation



## Corrective Action for Irritant Contact Dermatitis

- ① Consult Occupational Health
- ① Switch to gloves low in chemical residues (refer to data sheet for evidence of low irritant potential)
- ① Powder free
- ① Use caution around infectious agents
- ① Consult a dermatologist if symptoms persist



# Allergic Contact Dermatitis

- Also known as Type IV, Delayed Hypersensitivity or Chemical Allergy
- More than 2800 substances known to be contact sensitizers
- Allergic condition affecting genetically predisposed individuals
- Dose and rate dependent
- Glove-related allergic contact dermatitis can occur with latex and synthetic gloves
- Non glove-associated causes are fragrances (e.g. in shampoos, soaps etc), disinfectants and nickel
- Typically occurs 6-48 hours after allergen exposure
- More than 80% of diagnosed cases can be traced to accelerators



# Symptoms of Chronic Allergic Contact Dermatitis

- Dryness
- Scaling
- Peeling
- Pimples
- Reaction extends beyond the area of glove contact
- Cracking
- Inflammation
- Skin blisters





# Chronic form of Allergic Contact Dermatitis








## Corrective Action for Allergic Contact Dermatitis

- ① Consult Occupational Health
- ① Switch to gloves that are low in chemical allergens (refer to data sheet for evidence that gloves are low in chemical sensitizers)
- ① If specific chemical allergen has been identified, select a glove not containing this chemical



# Gloves -Allergen Exposure (Millions) for Allergic Contact Dermatitis

## 8 Changes/Day





# Natural Rubber Latex Allergy

- **Also known as also known as Immediate Type Hypersensitivity, Protein Allergy or Type I**
- **Allergic condition affecting genetically predisposed individuals**
- **Dose and rate dependent**
- **Can occur within one minute to an hour of exposure to allergen**
- **Potential for cross-reactivity with many common foodstuffs (e.g. apples, potatoes, tomatoes etc)**



# Symptoms of Natural Rubber Latex Allergy

- ① Urticaria
- ① Hives
- ① Asthma
- ① Anaphylaxis
- ① Hay fever symptoms



# Cross-Reactive Allergens

- Apples
- Avocados
- Bananas
- Celery
- Cherries
- Chestnuts
- Ficus
- Figs
- Grapes
- Kiwis
- Mangos
- Melons
- Passion Fruit
- Peaches
- Celery
- Pears
- Pistachios
- Potatoes
- Ragweed
- Strawberries
- Tomatoes
- Latex



## Natural Rubber Latex Allergy - Staff Must:

- Notify their supervisor and Occupational Health
- See an allergy specialist
- Wear synthetic gloves
- Work in a powder free environment
- Notify colleagues, dentist or GP
- Wear a medic alert bracelet
- Carry a source of adrenaline if prescribed





# TRIAGE FORMAT: GUIDE to DIFFERENTIATION of GLOVE REACTIONS

**EMPLOYEE COMPLAINS OF REACTION TO GLOVES**

Report to Occupational or Employee Health

This triage is only a guideline. Employer protocols regarding glove-associated reactions should be reviewed by your staff allergist and dermatologist.

Location: Hands or Arms Only

NO

Symptoms in addition to (or in lieu of) those on the hand/arm :

- Abdominal cramping, nausea
- Anaphylaxis
- Angioedema, pharyngeal swelling
- Blood pressure drop, tachycardia
- Diarrhea
- Dyspnea\*
- Headaches, disorientation
- Itching, burning eyes\*
- Respiratory distress, asthma\*
- Rhinitis\*

\* Note: If symptoms persist, but diagnostic test for immediate type hypersensitivity (Type I) to latex are negative, consider airborne irritants and chemical sensitizers (Type IV) of glove and non-glove origin

YES

Urticaria (Hives)

NO

How long have symptoms occurred?

YES

Recent (Acute)  
Symptom onset  
after glove removal

Long-term (Chronic)

Rapid

Itching, swelling,  
blisters, sores,  
cracking, redness  
or dry appearance

Delayed (6 to 72 hours)

Clustered bumps, itching (but  
upon scratching, painful),  
peeling, scales, blisters, red or  
dry appearance

Symptoms  
Appear beyond  
border of the glove

YES

Delayed onset

NO

NO

**IMMEDIATE TYPE HYPERSENSITIVITY  
(TYPE I)**

**IRRITATION**


**DELAYED TYPE HYPERSENSITIVITY  
(TYPE IV)**

**IRRITANT CONTACT  
DERMATITIS**





# Threshold for Clinical Symptoms – Natural Rubber Latex Allergy





# Natural Rubber Latex Products





# Employee OH review





# Glove barrier issues – Focus on vinyl

- Up to 60% in-use failure rate (Rego & Roley, 1999)
- Poor chemical resistance properties
- DEHP is a reprotoxin
  - Carcinogenic
  - Mutagenic
  - Reduces reproduction
  - PVC medical devices must be labelled “Containing phthalates”





# Glove Selection

<i>To reduce:</i>	<i>Choose gloves:</i>
Irritant Contact Dermatitis	<ul style="list-style-type: none"><li>➤ Low in chemicals</li><li>➤ Low in endotoxins</li><li>➤ Powder free</li></ul>
Allergic Contact Dermatitis	<ul style="list-style-type: none"><li>➤ Low in chemical contact sensitizers &amp; powder free</li></ul>
Natural Rubber Latex Allergy	<ul style="list-style-type: none"><li>➤ Low in protein</li><li>➤ Powder free</li></ul>
Occupational Asthma	<ul style="list-style-type: none"><li>➤ Powder free</li></ul>
In-use barrier failure	<ul style="list-style-type: none"><li>➤ Quality material appropriate to task</li></ul>






# Dermal Reactions

TYPE OF DERMAL REACTION	INCIDENCE OF OCCURRENCE
Irritant contact dermatitis	40%-60%* *Owenby, 1998
Allergic contact dermatitis	12%* *Gibbon, 2001
Natural rubber latex allergy	0.8%-7%* * Lebenbom-Mansour, 1997



## Disposable Glove Directives -which Directives are relevant?

- Workplace Directive 89/654/EEC
  - Employer obligations → Protect the employee
- PPE at Work Directive 89/656/EEC
  - Employers obliged to undertake risk assessment of hazards to employees
- Personal Protective Equipment Directive (PPE) 89/686/EEC
  - Manufacturer's obligations → Protect the wearer.
- Directive 93/42/EEC on Medical Devices (MDD).
  - Manufacturer obligations → Protect the patient.



THIS SUGGESTS THAT GLOVES FOR USE IN  
THE LABORATORY SHOULD BE  
« PPE » REGISTERED



# How to identify MDD Gloves?

- Underneath the CE mark, EN455 will usually feature providing easy identification.
- Non-sterile MDD Gloves are also often called Exam, Examination or Medical Gloves.
  - Highlighting their role in patient care.
- Non-sterile exam gloves are considered Class 1 Medical Devices.
  - This means gloves undergo a self certification process conducted by the manufacturer.
  - No independent validation of the test data by an external organization.

**CE**

**EN455**



## PPE Directive 89/686/EEC (cont) – summary of categorization

<b>PPE CATEGORY</b>	<b>RISK</b>	<b>CERTIFICATION / INSPECTION</b>	<b>CE MARKING</b>
PPE Category I: Simple Design	Minimal	Self-certification	<b>CE</b>
PPE Category II. Intermediate Design	Intermediate	Notified Body (product testing)	<b>CE</b>
PPE Category III. Complex Design	Irreversible or mortal	Notified Body (product testing+ quality audit)	<b>CE</b> <b>0321</b>



# INTERPRETATION OF EU REGULATIONS

## PROTECTION AGAINST VIRUSES

- EN374-2:2003 is the reference standard
  - Based on the water leak or air leak procedure to demonstrate that gloves offer an effective barrier to micro-biological hazards
  - NB: EN374-1:2003 Article 3.2 Explains that the micro-organism resistance definition extends only to bacteria and fungi. The test does NOT apply to protection against viruses.

THE VIRAL PENETRATION TEST  
(ASTM F1671)  
COULD BE THE RIGHT SOLUTION





# SHIELDskin™ ORANGE NITRILE™ Gloves


## WATER LEAK ASSAY

- **Impact:** Leaks and defects on gloves impact *personal safety* and *process integrity* from human-borne contamination
- **Test Method:** **EN374-2:2003** “Protective Gloves Against Chemicals and Micro-organisms – Part 2: Determination of Resistance to Penetration” is suitable for demonstrating that gloves form an effective barrier to bacteria and fungi but *does not apply to viruses*
- **Limits:** To be considered micro-organism resistant gloves must achieve at least an AQL of <1.5. All SHIELDskin™ gloves must achieve an AQL of <0.65 *making them more than twice as effective*
- **Testing:** is conducted for every manufactured batch of gloves. Additionally regular external auditing is carried out by a *Notified Body* (indicated by the 4 digits under the CE mark) as part of the PPE Cat III certification.



## VIRAL PENETRATION TEST

- **Impact:** Users handling live viruses need to know if their glove can protect them from accidental spills as well as provide a suitable barrier to protect process integrity. This test evaluates the quality of the glove film that cannot be detected using Water Leak Assay
- **Test Method:** **ASTM F1671-07** “Test Method for Resistance of Materials Used in Protective Clothing to Penetration by Blood-Borne Pathogens” using the non-hazardous Phi-X174 Bacteriophage as a penetration test system
- **Limits:** No viruses should go through
- **Industry standard method:** Yes
- **Testing:** Required for every new glove or glove formulation, and tests are done on as needed basis. Performed at an outside laboratory such as Nelson Laboratories





# Choosing the right glove

- 1) What is the principal intended purpose:
  - o Personal Protection
  - o Process Protection
  - o Patient Protection
- 2) What is the risk
  - o Minimal Risk => PPE Cat I i.e. easily reversible
  - o Irreversible / mortal risk => PPE CATEGORY III
- 3) How do I protect my hand from the wearing of gloves
  - o Glove associated irritation
    - o Powder free
    - o Low in chemical residues (viz. primary skin irritation test)
  - o Glove associated Type IV reactions
    - o Powder free
    - o Low in chemical contact sensitizers (viz. Buehler sensitization test, Modified Draize test & specific chemical sensitizer test)
  - o Glove associated Type I reactions
    - o latex-free gloves
- 4) Which material do I select – vinyl, nitrile or latex?
  - o Long cuff Nitrile preferred
  - o Ensure proper fit

